


MINI
S
C
E
N
E
S

SHENYUN

Sterile Single-Use Abdominal Wall Suturing Device

An innovative medical device integrating
puncture and suturing functions,
specifically designed for minimally invasive surgery.



PRODUCT COMPOSITION

Sterile Single-Use Abdominal Wall Suturing Device

The disposable abdominal wall suture kit consists of a puncture device (C) and an abdominal wall suture device (F). Pneumoperitoneum needle (Q) and specimen collection bag (D) are optional accessories.

- Puncture device: The sizes are available in eight specifications: 5, 10, 12, 15, 5A, 10A, 12A, and 15A. The letter "A" indicates the long model.

The absence of letters indicates regular styles.

- Pneumoperitoneum needle: It is available in two specifications, 120 and 150, based on the length of the cannula.
- Abdominal wall suture device: It is divided into four specifications according to the different diameters of the needle tubes, namely 1.4, 1.6, 1.8 and 2.0.
- Specimen collection bag: Maximum capacity 250mL.



CLINICAL APPLICATION



General Surgery

Gastrointestinal tract, liver, gallbladder and pancreas, colorectal tract, urinary tract, etc.



Thoracic Surgery

Thoracic surgeries (such as lung resection, esophageal surgery), etc.



Urology

Nephrectomy, partial nephrectomy, total cystectomy, prostate surgery, etc.



Obstetrics and Gynecology

Cesarean section, hysterectomy, myomectomy, ovarian tumor resection, etc.



Bariatric surgery

Single-department products are recommended as the first choice for in-hospital departments.

Suzhou Shenyun Medical Equipment Co., Ltd.
Add: No. 666, Yinzang Road, Linhu Town, Suzhou City 215104, China
P: + 86 182 2153 2164
E: Jolin@symedical.com.cn
www.shenyunmedical.com



Scan to watch the steps for using the device

■ PRODUCT ADVANTAGES

Advantage

Integrating a puncture device and a fascia suture device into one function, it saves the number of instruments used.

Simple operation

Threading the needle and threading the thread is simple and convenient to operate, saving operation time



High efficiency

Precise positioning enables faster capture of stitches, allowing for blind sewing.

Prevent puncture hole hernia

Full-thickness suture leaves no dead space, effectively preventing puncture hole hernia after laparoscopic surgery.

■ SPECIFICATION AND SIZE

Abdominal Wall Suturing Device

| Specification | Needle length (mm) | Needle diameter (mm) |
|---------------|--------------------|----------------------|
| 2.0 | 160 | 2.1 |
| 1.8 | 160 | 1.8 |
| 1.6 | 160 | 1.6 |
| 1.4 | 160 | 1.4 |

Veress Needle

| Specification | Inner diameter (mm) | Total Length (mm) |
|---------------|---------------------|-------------------|
| 120 | 1.7 | 120 |
| 150 | 1.7 | 150 |

Specimen Retrieval Bag

| Specification | Bag opening diameter (mm) | Maximum capacity (ml) |
|---------------|---------------------------|-----------------------|
| 1.0 | 95 | 250 |

Trocar

| Specification | Length (mm) | Inner diameter (mm) |
|---------------|-------------|---------------------|
| 5 | 110 | 5.8 |
| 10 | 110 | 10.9 |
| 12 | 110 | 12.9 |
| 15 | 110 | 15.6 |
| 5A | 150 | 5.8 |
| 10A | 150 | 10.8 |
| 12A | 150 | 12.5 |
| 15A | 150 | 15.6 |